

Optimal Results Require Proper Specimen Collection



Collection Videos and Guides for Clinicians



**Easily accessible
video resources**



**Help reduce specimen
collection errors**



**Feel comfortable with
the collection process**

View step-by-step specimen collection guides and videos for optimal results.

ThinPrep® Vial

Out-of-the-vial Testing:

- ThinPrep® Pap test
- Aptima® HPV
- Aptima® HPV 16 18/45 Genotyping
- Aptima® Combo 2® assay for CT/NG
- Aptima® Trichomonas vaginalis assay

Clinician Collection Site:

- Cervix



Aptima® Multitest Swab

Assays:

- Aptima® Combo 2 (CT/NG)*
- Aptima® Mycoplasma genitalium*
- Aptima® Trichomonas vaginalis*
- Aptima® BV*
- Aptima® CV/TV*
- Aptima® HSV 1 & 2

 Patient or self-collection devices are approved for all Aptima assays, except Aptima HSV 1&2.

Clinician Collection Site:

- Vaginal*
- Penile meatal*
- Throat
- Rectal
- Genital lesion

 Patient or self-collection available for vaginal or penile meatal sites.



Aptima® Urine

Assays:

- Aptima® Combo 2 (CT/NG)
- Aptima® Mycoplasma genitalium
- Aptima® Trichomonas vaginalis



Aptima® Unisex Swab

Assays:

- Aptima® Combo 2 (CT/NG)
- Aptima® Mycoplasma genitalium
- Aptima® Trichomonas vaginalis

Clinician Collection Site:

- Endocervical
- Urethral



Refer to individual assay package inserts for cleared specimen types and performance claims.^{1,9}

* Patient collection available for this device or source.

References: 1. ThinPrep Pap Test Quick Reference Guide.#DS-05867. Marlborough, MA: Hologic, Inc.; 2017. 2. Aptima HPV Assay. US package insert AW-12820, Rev 004. San Diego, CA: Hologic, Inc.; 2020. 3. Aptima HPV 16 18/45 [package insert]. AW-26281, Rev 001. San Diego, CA: Hologic, Inc.; 2023. 4. Aptima Combo 2 Assay. Pather. US package insert AW-25929-001. Hologic, Inc.; 2023. 5. Aptima Trichomonas vaginalis Assay. Pather. US package insert AW-14170-001. Hologic, Inc.; 2023. 6. Aptima Mycoplasma genitalium Assay. Pather. US package insert AW-17946-001. Hologic, Inc.; 2023. 7. Aptima BV Assay. Pather. US package insert AW-23712-001. Hologic, Inc.; 2023. 8. Aptima CV/TV Assay. Pather. US package insert AW-23713-001. Hologic, Inc.; 2023. 9. Aptima Herpes Simplex Viruses 1 & 2 Assay. Pather. US package insert AW-15636-001. Hologic, Inc.; 2021.

MISC-07761-001 Rev. 003 © 2023 Hologic, Inc. All rights reserved. Hologic, ThinPrep, Aptima and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to diagnostic.solutions@hologic.com.

Collection and Video Guides



Scan to discover more

